

• Kit Specifications:

	Reagent/Quantity	Storage
Cat. No.: CK0011	R.1: 2x 50 ml R.2 1x20ml total 120 ml	2-8°C
Cat. No.: CK0111	R.1: 4x 50 ml R.2 1x40ml total 240 ml	2-8°C
Cat. No.: CK0017	R.1: 2 x 50 ml R.2 1x 20ml total 120 ml	2-8°C
Cat. No.: CK0117	R.1: 4x 50 ml R.2 1x40ml total 240 ml	2-8°C

• Intended Use:

In Vitro Diagnostic reagent pack for the quantitative determination of Creatine Kinase-NAC in human serum/plasma on automated and semi-automated photometric systems.


• Summary and Explanation:

The Creatine Kinase is the dimer enzyme which exist in 4 types. The one iso enzyme is in mitochondria and three iso enzymes in cytoplasm: CK-MM (muscular enzyme), CK-BB (brain enzyme) and CK-MB (heart enzyme). The CK enzyme activity contain the all iso enzymes and Total CK. The healthy person CK activity is more related to CK-MM and the other iso enzymes activity limited. The CK enzyme which increase in heart and skeleton cell injuries so the measurements of this enzyme is critical. The measurement of CK-MB activity joint with CK-NAC value is an exclusive test for diagnostic of heart muscle damage and determine the early (previous occur) diagnostic of myocardial infarction. The CK-NAC will increase activity 3 to 8 hours before infarction. Also, the enzyme activity related to others disease such as skeleton muscle injury, heart, CNS and thyroid. Moreover, the increase activity of this enzyme is related to brain vascular disease, acute muscle dystrophy, dermatomyositis and electric shock.

Principle of the Method:

In this method the creatine kinase catalyze creatine phosphate to creatine and ATP. Furthermore, hexokinase change the ATP and glucose to G6P and ADP. In the final step the G6P-DH convert G6P with NADP+ to 6-phosphogluconate and NADPH + H+. The speed formation of NADPH is measurement on 340 nm which is related to amount of CK-NAC in the sample.

CK (pH6.7)



• Reagent Preparation and Stability:

Reagent is ready for use.

Before use, mix reagent by gently inverting each bottle.

Reagent is stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Do not freeze and protect from light. For assay the mixture should mix by with 1 volume of R2 and 5 volumes of R1. R1 is transparent and colorless also R2 is same as well. The mixture ready to use, stability at 2-8°C is 15 days. Avoid direct sunlight.

Waste Management:

Refer to local legal requirements for chemical disposal regulations.

Warning: Handle waste as potentially biohazardous material.

Dispose of waste according to accepted laboratory instructions and procedures.

• Warnings and Precautions:

For In Vitro Diagnostics Use Only.

For Professional Use Only.

Carefully read instructions for use.

In case of serious damage to the bottle or cap, resulting in product leakage or contamination, do not use the reagent pack and contact your distributor.

Take all necessary precautions required when handling laboratory reagents.

Do not use components past the expiry date stated on the Bottles.

Do not interchange caps among components as contamination may occur and compromise test results. For diagnostic purposes, the results should always be

assessed with the patient's medical history, clinical examinations and other findings.

Any serious incident related to the product must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is located.

• Type of Specimen:

Use fresh serum non-heamolysied. There are some references recommended to use plasma for this assay as the speed of change in plasma is really fast and unpredictable.

The serum/plasma should have collected immediately after blood collection. Avoid any contamination. The CK-NAC stability in the sample of serum/ plasma at 2-8°C remain till 7 days but in at -20°C stable for 30 days. Avoid any contamination and direct sunlight. The sample collector door should be close tightly.

Required but not Supplied:

General chemistry calibrator from TKS or other valid calibrators.

General chemistry control Level 1 & 2 from TKS or other valid controls.

Saline solution 0.9 % NaCl

General laboratory equipments

• Notes:

Carefully read instructions for use.

It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with hydrochloric acid 1 N and then thoroughly rinsed it with distilled water.

Use clean disposable pipette tips for its dispensation.

Disposable and glassware material used must be free of metals, ions and detergents.

Performance Characteristics:

Performance results can vary with the instrument used.

Data obtained in each individual laboratory may differ from these values.

Maximum determination in this assay is 1800U/L

LOD: 5 U/L

For samples with a higher concentration (1800U/L), dilute 1:1 with 0.9 % NaCl and re-assay. Multiply result by 2.

Precision:

Intra Assay-Within run CK-NAC

Sample	n	Mean (U/L)	SD (U/L)	CV (%)
1	20	169	1.50	0.89
2	20	383	1.29	0.34

Inter Assay-Between run CK-NAC

Sample	n	Mean (U/L)	SD (U/L)	CV (%)
1	20	177	2.5	1.41
2	20	408	3.43	0.84

• Accuracy:

Results obtained using BIOMEDIC reagents (y) did not show systematic differences when compared with other commercial reagents (x).

Correlation coefficient (r): 0.986

Regression equation: Y = 0.997 (X) + 5.765 U/L

The results of the performance characteristics depend on the analyzer used.

• Interfering Substances:

the less concentration of below items are not interfering in this assay.

Bilirubin (mixed isomer)	Less than 10% interference up to 400 μmol/L Bilirubin
Lipaemia	Less than 10% interference up to 5 g/L Intralipid
Haemolysis	Less than 10% interference up 1.25g/L Hemoglobin.

• Reference Values:



CK men	>80U/L at 25°C	>130U/L at 30°C	>190U/L at 37°C
CK women	>70U/L at 25°C	>110U/L at 30°C	>167U/L at 37°C
CK-MB	>10U/L at 25°C	>15U/L at 30°C	>24U/L at 37°C
CK-MB activity amounts for 6-25% of the total CK activity.			

If the three situations occurred in the same time for patient, the high possibility is heart attacked proceed for patient.

CK at 37°C	
Women	24-167 U/L
Men	24-190 U/L
New born	468-1200 U/L
New born till 5 days	195-700 U/L
New born till 6 months	41-330 U/L
New born more than 6 months	44-229 U/L

Each laboratory should establish its own expected values. The creatine kinase total results should always be reviewed with the patient's medical examination and history.

Assay Procedure:

Allow reagents to reach working temperature before using.

A proportional variation of the reaction volumes indicated does not change the result

Assay conditions:

334-365 (340)nm	Wavelengths
30°C (25°C or 37°C)	Incubation Temperature
1 cm	Cuvette

Adjust the instrument to zero with distilled water.

Control/Sample/Calibrator	Blank	
1000 µl	1000 µl	R1
40 µl	-	Control/Sample/Calibrator
Gently mix and incubate for 5 minute at 37°C. Then added R2.		

200 µl	200 µl	R2
Gently mix and incubate for 2 minute at 37°C. Measure the sample and calibrator. Turn on the counter exactly 1,2 and 3 min after start reaction measure the absorbance. ΔOD/min.		

For single reference following the below table

Control/Sample/Calibrator	Blank	
1000 µl	1000 µl	R
40 µl	-	Control/Sample/Calibrator
Gently mix and incubate for 2 minute at 30°C (25°C or 37°C). Measure the sample and calibrator. Turn on the counter exactly 1,2 and 3 min after start reaction measure the absorbance. ΔOD/min.		

• $\Delta A/\text{min} = \Delta A/\text{sample or calibrator} - \Delta A/\text{min Blank}$

Calculations:

$$\text{CK-NAC (U/L)} = \frac{\text{OD Sample} \times \text{Cal Conc. (U/L)}}{\text{OD calibrator}}$$

Conversion units:

$$\text{CK-NAC (}\mu\text{Kat/l)} = \text{CK-NAC (U/L)} \times 0.0167$$















TOSE'E KIMIA SA'ADAT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

References:

- Stein W. Med Weit. 1985, 36:572 & Burtis CA, Ashwood ER. Tietz Fund. Of Clin. Chem. 5th ed. 30-54, 352-390 and 974-975.
- GuderWG, NarayananS, WisserH, ZawtaB. List of Anal; Preanal Variables. From the Patient to the Laboratory. Darmstadt: GIT Verlag 1996.
- Wurzberg U, Hennrich N, Lang H, Prellwitz W, Neumeier D, Knedel M. Klin. Wschr. 1976: 54 and 357.
- Szasz G, Busch EW. Third European Congress of Clinical Chemistry, Brighton, England, 3-8 June 1979 (abstract).

Symbols:

The following symbols are used in the labelling of TOSE'E KIMIA SA'ADAT systems:

 In Vitro Diagnostics	 Contains sufficient for <n> tests
 Batch Code	 Temperature limit
 Catalogue No.	 Consult instruction for use
 Expiry Date	 Caution
 Date of Manufacture	 Keep dry
 Manufactured by	 This way up
 Biological Risks	 Keep away from sunlight

Rev 01: Issued on 20 February 2023

